

Accessories

High-Pressure Injector Lines

HP INJECTOR LINES 1200 PSI - NON-BRAIDED

- Non-braided
- Ideal for all power injector applications, offered in a variety of lengths
- Available with or without High-Pressure Rotating Adaptor
- Contents per case: 25

| | LENGTH (cm) | DESCRIPTION | CATALOGUE NUMBER |
|--|----------------|----------------------|---------------------|
| | 25 | with rotator | 12-000813 |
| | 50 | with rotator | 12-001771 |
| | 75 | with luer-lock, male | 12-000237 |
| | 75 | with rotator | 12-001733 |
| | 120 | with luer-lock, male | 12-000608 |
| | 120 | with rotator | 12-001734 |

HP INJECTOR LINES 1200 PSI - BRAIDED

- Braided
- Combines the high pressure capacity with superior flexibility
- Available with or without High-Pressure Rotating Adaptor
- Contents per case: 25

| | LENGTH (cm) | DESCRIPTION | CATALOGUE NUMBER |
|--|----------------|----------------------|---------------------|
| | 25 | with luer-lock, male | 12-001743 |
| | 50 | with rotator | 12-001738 |
| | 75 | with luer-lock, male | 12-001741 |
| | 75 | with rotator | 12-001737 |
| | 120 | with luer-lock, male | 12-001740 |
| | 120 | with rotator | 12-001736 |

Manifolds 1050 PSI

- Continuous large internal diameter (2.3 mm diameter) ideal for special procedures requiring small catheters.

| | DESCRIPTION | CATALOGUE NUMBER |
|--|------------------------------------|---------------------|
| | 2-port, with rotator, right-handed | 12-001759 |
| | 3-port, with rotator, right-handed | 12-001765 |
| | 3-port, with rotator, right-handed | 12-001767 |



Accessories

High-Pressure Injector Lines

Manifolds 550 PSI

| | DESCRIPTION | CATALOGUE NUMBER |
|--|------------------------------------|------------------|
| | 2-port, with rotator, right-handed | 12-001727 |

Pressure Monitoring Lines

- Clear and high durometer for accurate measurements and optimum dynamic response

| | LENGTH (cm) | DESCRIPTION | CATALOGUE NUMBER |
|--|-------------|--------------|------------------|
| | 15 | male/female | 12-000235 |
| | 30 | male/female | 12-000836 |
| | 90 | male/female | 12-001776 |
| | 120 | male /female | 12-000236 |
| | 150 | male/female | 12-000210 |



Accessories

Extensions / Connectors / Inflation Systems

CATHETER EXTENSIONS

- Flexible, translucent polyurethane
- Braided
- Designed for 1000 psi
- Luer lock hubs
- Male-female
- 5 units per package



| | LENGTH (cm) | INNER DIAMETER (mm) | CATALOGUE NUMBER |
|--|----------------|------------------------|---------------------|
| | 45 | 1.8 | 502-100D |
| | 65 | 1.8 | 502-101D |
| | 100 | 1.8 | 502-102D |

VALVES & HAEMOSTATIC Y-CONNECTORS

- Available for different wire diameters
- Smooth turning Tuohy-Borst valve
- With rotating adaptor
- Available as a PTCA valve accessory kit with either a tool or a puncture device and a torque device

| | CASE CONTENT | DESCRIPTION | CATALOGUE NUMBER |
|--|--------------|---|---------------------|
| | 25 | Large Bore haemostatic Y-connector for catheter up to 9F (Torque device - Guidewire Insertion Tool - Haemostatic Y-connector for Catheters up to 9F) | 12-001915 |
| | 25 | Haemostatic Y-connector for catheters up to 9F | 12-001855 |
| | 25 | Torque device | 12-001624 |
| | 25 | Twin interchange Y-connector for catheters up to 9F | 12-001856 |

INFLATION SYSTEMS

| | CASE CONTENT | DESCRIPTION | CATALOGUE NUMBER |
|--|--------------|---|---------------------|
| | 6 | Disposable Inflation Device 30ml pressure control syringe with manometer measuring pressures between 0 and 30 atm (incl 3-way stopcock) | 12-003704 |



Accessories

Control Syringes

CONTROL SYRINGES

- Smooth operation
- Clear, glass-like plastic syringe for clarity
- Contents per case: 25



| | VOLUME (ml) | DESCRIPTION | STOPPER (ml) | CATALOGUE NUMBER |
|--|----------------|----------------------|-----------------|---------------------|
| | 10 | with luer-lock, male | 0.5 | 12-000249 |
| | 10 | with rotator, male | 0.5 | 12-000304 |



Accessories

Angiographic Needles

ONE-PIECE NEEDLES

- Thin wall cannula
- Luer hub
- 10 units per package



| CANNULA OD | DESCRIPTION | MAXIMUM GUIDEWIRE (") | LENGTH (cm) | CATALOGUE NUMBER |
|------------|-------------|--------------------------|----------------|---------------------|
| 18G | Thin wall | 0.038 | 7 | 12-004943 |
| 19G | Thin wall | 0.032 | 3.8 | 12-004925 |
| 21G | Thin wall | 0.021 | 3.8 | 12-004924 |

ONE-PIECE NEEDLES

- Thin wall cannula
- Luer hub
- 25 units per package



| CANNULA OD | DESCRIPTION | MAXIMUM GUIDEWIRE (") | LENGTH (cm) | CATALOGUE NUMBER |
|------------|------------------------------|--------------------------|----------------|---------------------|
| 18G | Thin wall with snap-on wings | 0.038 | 7.0 | 502-652 |
| 19G | Thin wall with snap-on wings | 0.032 | 7.0 | 502-654* |
| 20G | Thin wall with snap-on wings | 0.021 | 5.0 | 502-655* |
| 20G | Thin wall with snap-on wings | 0.025 | 5.0 | 502-659* |

TWO-PIECE NEEDLES

- Thin wall cannula
- Luer hub

- Solid stylet
- 25 units per package



| CANNULA OD | DESCRIPTION | MAXIMUM GUIDEWIRE (") | LENGTH (cm) | CATALOGUE NUMBER |
|------------|------------------------------|--------------------------|----------------|---------------------|
| 18G | Thin wall with snap-on wings | 0.038 | 7.0 | 502-657 |
| 19G | Thin wall with snap-on wings | 0.035 | 7.0 | 502-658 |

Accessories

Manifolds

250 PSI MANIFOLDS

| | | DESCRIPTION | CATALOGUE NUMBER |
|--|-----|------------------------------------|------------------|
| | ON | 2-port, with rotator, right-handed | 12-004948 |
| | OFF | 2-port, with rotator, right-handed | 12-004951 |
| | OFF | 3-port, with rotator, right-handed | 12-004952 |