


SLEEK™ OTW

.014" Over-the-Wire PTA Dilation


Broad portfolio of sizes including 280 mm length

Excellent Pushability and a Small Crossing Profile.


- Choose from 1.25-5 mm diameters and 15-280 mm lengths
- Long 150 cm working catheter lengths



Excellent Pushability for More Distal Revascularization¹



Small Crossing Profile for Tighter Lesions¹



Product Benefits:

- **Excellent pushability** with a **small crossing profile** to facilitate the treatment of the most challenging and distal infrapopliteal lesions
- Offers a **broad portfolio** enabling you to treat patients with critical limb ischemia more easily and with fewer inflations

1. Cordis Data on file.



Ordering Information

SLEEK™ Over-the-Wire 150 cm Delivery System

Catalog Number	Balloon Diameter, mm	Balloon Length, mm	Rated Burst Pressure, atm	Sheath Fit, F
426-1201X	1.25	15	14	4
426-1504X	1.5	40	16	4
426-1508X	1.5	80	16	4
426-1510X	1.5	100	16	4
426-2002X	2	20	16	4
426-2004X	2	40	16	4
426-2008X	2	80	15	4
426-2010X	2	100	15	4
426-2012X	2	120	15	4
426-2015X	2	150	15	4
426-2022X	2	220	15	4
426-2028X	2	280	15	4
426-2502X	2.5	20	16	4
426-2504X	2.5	40	16	4
426-2508X	2.5	80	15	4
426-2510X	2.5	100	15	4
426-2512X	2.5	120	15	4
426-2515X	2.5	150	15	4
426-2522X	2.5	220	15	4
426-2528X	2.5	280	15	4
426-3002X	3	20	16	4
426-3004X	3	40	16	4
426-3008X	3	80	15	4
426-3010X	3	100	15	4
426-3012X	3	120	15	4
426-3015X	3	150	15	4
426-3022X	3	220	15	4
426-3028X	3	280	15	4
426-3502X	3.5	20	16	4
426-3504X	3.5	40	16	4
426-3508X	3.5	80	15	4
426-3510X	3.5	100	15	4
426-3515X	3.5	150	15	4
426-3522X	3.5	220	15	4
426-3528X	3.5	280	15	4
426-4002X	4	20	16	4
426-4004X	4	40	16	4
426-4008X	4	80	15	4
426-4010X	4	100	15	4
426-4015X	4	150	15	4
426-4022X	4	220	15	4
426-4028X	4	280	15	4
426-5002X	5	20	14	4
426-5004X	5	40	14	4
426-5008X	5	80	13	4
426-5010X	5	100	13	4
426-5015X	5	150	10	4
426-5022X	5	220	10	4
426-5028X	5	280	10	4

For Healthcare Professionals Only.

Important Information: Prior to use, refer to the "Instructions for Use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings, and precautions. As part of the Cordis policy of continuous product development, we reserve the right to change product specifications without prior notification. For healthcare professionals only.

© Cardinal Health 2017. All right reserved. CORDIS, the Cordis LOGO, and SLEEK are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. EU2352 06/17.

Cordis
A Cardinal Health company